




basic kit

product information



Flavin-mononucleotide-based Fluorescent Protein (FbFP)

evoglow[®] basic kit

Catalog No. evo-2.1.010

Quantity 20 µg each

General Information

Fluorescent reporter proteins are valuable noninvasive molecular tools for *in vivo* real-time imaging of living cells and tissues as well as *in vitro* fluorescence labeling. A major drawback of existing GFP-like reporter proteins is their strict requirement for molecular oxygen as a cofactor for the synthesis of their respective chromophores. Therefore, the application in anaerobic systems is not possible (see Fig. 1).

The flavin mononucleotide (FMN)-based fluorescent proteins (FbFPs) from the evoglow[®] series were developed to overcome these restrictions. The evoglow[®] proteins can be used as fluorescent reporters in both aerobic and anaerobic biological systems.

Kit content: The evoglow[®] basic kit contains the plasmids pGLOW-Bs1, pGLOW-Bs1-stop, pGLOW-Bs2, pGLOW-Bs2-stop, pGLOW-Pp1 and pGLOW-Pp1-stop. The genes are cloned *SacI/PstI* or *SacI/BamHI* into the multiple cloning site of pUC18 for individual subcloning. Each gene is available in a fusion-version - without stop codon - and in an expression-version containing the TGA stop codon.

Origin: The evoglow-Bs1 and Bs2 genes were cloned from the Gram positive bacterium *Bacillus subtilis*. The Bs1-sequence is the wild-type codon-usage, whereas Bs2-sequence has been adjusted to *Escherichia coli* codon usage and is truncated, i.e. comprises the LOV-domain only. The evoglow-Pp1 gene was cloned from the Gram negative bacterium *Pseudomonas putida*.

Formulation: pGLOW plasmids are delivered as 20 µg plasmid DNA dissolved in regular 10 mM Tris/HCl-buffer, pH 8,5.

Storage: pGLOW plasmids can be stored at -20°C.

Intended use: For research use only. Not for clinical diagnosis.

References: Drepper, T. *et al.* (2007) *Nat. Biotechnol.* 25: 443-445.





Fig. 1: *in vivo*-fluorescence of evoglow[®]-proteins compared to that of conventional Yellow Fluorescent Protein (YFP).



evoglow®- Anaerobic Fluorescent Proteins

The designation evoglow® comprises a novel type of fluorescent proteins containing a Flavin-mono-nucleotide-based cofactor. They are thus capable of developing bright cyan-green fluorescence even in the complete absence of oxygen. They fold rapidly *in vivo* and thus enable immediate detection of expression to study anaerobic inter- and intracellular processes. Since

their fluorescence intensity as well as the excitation and emission spectra of the evoglow® proteins lie within the same range as those of other fluorescent proteins, the same instrumental setup can be used for their detection and analysis as for other fluorescent proteins. Recommendations for instruments and filtersets are given in section *Instrumental Setup*.

Physical Parameters:

Characteristics	evoglow®-Bs1	evoglow®-Bs2	evoglow®-Pp1
Fluorescence color	cyan-green	cyan-green	cyan-green
Excitation max. (nm)	450	450	450
Emmission max. (nm)	495	495	495
Quantum yield	0,4	0,39	0,17
Extinction coefficient (M ⁻¹ cm ⁻¹)	12,500	12,500	12,500
Brightness*	5	4.9	2.1
pKa	n.d.	n.d.	n.d.
structure	monomer	dimer	dimer
Photostability	+	++	+++
Molecular weight	ca. 33kDa	ca. 19kDa	ca. 19kDa
Amino acids	262	137	148

Instrumental Setup Recommendations


evoglow® proteins can be excellently detected by using the following instrumental setup:

- microscope: Zeiss Axioplan 1.
- object lens for fluorescence imaging: Plan APOCHROMAT 100x Öl DIC mit numerical aperture 1,4 (object lens "Plan-Apochromat" 100x/1,40 Oil DIC) (Immersion oil 518F)
- filterset 37 (Ex 450/50, BS 480, Em 510/50)
- lamp: HBO50
- camera: Zeiss AxioCam MRm (12 Bit, monochrome)
- standard imaging software: AxioVision Rel. 4.6.
- for light microscopy: Zeiss Plan NEOFLUAR 100x Öl phase contrast objective (Phase 3) numerical aperture 1,3 (object lens "Plan-Neofluar" 100x/1,30 Oil Ph3)

pGLOW-Bs1 and pGLOW-Bs1-stop

Plasmid Information

Plasmid Map



In pGLOW-Bs1-stop the evoglow-Bs1-gene includes a termination codon indicated in red

```

5' >> tac gaa ttc gag ctc cat atg... evoglow-Bs1-stop ...atg tga gga tcc tct aga gtc gac ctg cag gca tg...>> 3'
   >> atg ctt aag ctc gag gta tac... ..tac act cct agg aga tct cag ctg gac gtc cgt ac...>>
  
```


Limitation of use:

The plasmids distributed by company evocatal GmbH under the trade mark "evoglow®" are subject to intellectual property right applications. Company evocatal GmbH commits the plasmids and genes to the customer solely for the purpose of own scientific, non-commercial investigation or own, non-commercial research and development purposes. Any commercial use requires the conclusion of a separate license contract. For commercial customers (i.e. companies or any for-profit-organisations) the permission of use for the above said research purposes of the evoglow genes and plasmids is limited to six (6) months, starting from the date of purchase. For detailed license information please contact evocatal GmbH, Germany. Furthermore, the modification of the evoglow-gene-sequences and of the plasmids is assertive prohibited. By the purchase of the plasmids without notice of defect within a term of 14 days beginning with the day of delivery, the customer agrees to the limitation of use above.

Exclusion of liability:

Company evocatal GmbH commits the plasmids distributed under the trade mark "evoglow®" to the customer solely for use in microorganism. The plasmids are not intended for the use in human medical or veterinarian medical diagnostic purposes, or in the animal or human body; any liability of the company evocatal GmbH is excluded. The committed plasmids are classified to security level S1. The national safety regulations for handling genetic engineered microorganism have to be considered and fulfilled when handling the plasmids. By the purchase of the plasmids without notice of defect within a term of 14 days beginning with the day of delivery, the customer agrees to the above exclusion of liability and indemnifies company evocatal GmbH from any liability for the use of the plasmids beyond the intended application.

Fluorescence Spectra:



Technical Data:

excitation max.:	449 nm
emission max.:	495 nm (λ_{ex} : 450 nm)
molar extinction coefficient (FMN):	$\epsilon = 12,500 \pm 500 \text{ M}^{-1} \text{ cm}^{-1}$ (at 450 nm)
quantum yield:	$Q_f = 0,38$

Sequence Information:

evoglow-Bs1 nucleotide sequence (786 bp)

```

1  ATGGCTAGTT  TTCAATCATT  TGGGATACCA
31  GGACAGCTGG  AAGTCATCAA  AAAAGCACTT
61  GATCACGTGC  GAGTCGGTGT  GGTAATTACA
91  GATCCGCAC  TTGAAGATAA  TCCTATTGTC
121  TACGTAAATC  AAGGCTTTGT  TCAAATGACC
151  GGCTACGAGA  CCGAGGAAAT  TTTAGGAAAG
181  AACGCACGCT  TCTTACAGGG  GAAACACACA
211  GATCCTGCAG  AAGTGGACAA  CATCAGAACC
241  GCTTTACAAA  ATAAAGAACC  GGTCACCGTT
271  CAGATCCAAA  ACTACAAAAA  AGACGGAACG
301  ATGTTCTGGA  ATGAATTAAA  TATTGATCCA
331  ATGGAAATAG  AGGATAAAAC  GTATTTTGTC
361  GGAATTCAGA  ATGATATCAC  CAAGCAAAAA
391  GAATATGAAA  AGCTTCTCGA  GGATTCCTTC
421  ACGGAAATTA  CTGCACTTTC  AACTCCTATT
451  GTCCCGATTC  GCAATGGCAT  TTCGGCTCTT
481  CCGCTAGTCG  GAAACCTGAC  AGAGGAGCGA
511  TTTAATCCA  TCGTTTGAC  ATTGACGAAT
541  ATCTTATCAA  CATCCAAAGA  TGATTATTTG
571  ATCATTGATT  TATCCGGATT  GGCCCAAGTG
601  AACGAACAAA  CGGCCGACCA  AATTTTCAAG
631  CTGAGCCATT  TGCTGAAATT  GACCGGAACT
661  GAGTTAATCA  TTACTGGCAT  TAAGCCTGAA
691  TTGGCTATGA  AAATGAATAA  ACTGGATGCC
721  AATTTTTCGT  CGCTGAAAAC  ATATTCAAAT
751  GTAAAGGATG  CCGTTAAAGT  GCTTCCGATT
781  ATG  TGA term. codon in pGLOW-Bs1-stop only

```

evoglow-Bs1 amino acid sequence (261 aa)

```

1  MASFQSGFIP  GQLEVIKKAL  DHVRVGVVIT
31  DPALEDNPIV  YVNQGFVQMT  GYETEEILGK
61  NARFLQ GKHT  DPAEVDNIRT  ALQNKEPVTV
91  QIQNYKKDGT  MFWNELNIDP  MEIEDKTYFV
121  GIQNDITKQK  EYEKLLDSL  TEITALSTPI
151  VPIRNGISAL  PLVGNLTEER  FNSIVCTLTN
181  ILSTSKDDYL  IIDLSGLAQV  NEQTADQIFK
211  LSHLLKLTGT  ELIITGIKPE  LAMKMNKLDA
241  NFSSLKTYSN  VKDAVKVLP  I  M

```

Note:

The nucleotide sequence of **evoglow-Bs1-stop** is identical to that of **evoglow-Bs1** but includes a termination-codon at the 3'-end as indicated in the adjacent sequence. The plasmid base pair numbering is hence shifted 3 positions in the 'stop'-constructs.


All plasmids in this kit are derived from pUC18 and replicate only in *E. coli* and closely related microorganisms.

Plasmids are not suitable for expression.

pGLOW-Bs2 and pGLOW-Bs2-stop

Plasmid Information

Plasmid Map



restriction analysis information:

enzyme	fragment size (bp)
Alw441	497
	1246
	1331

further restriction sites at (bp)

EcoRI	851
HindIII	399
NdeI	183 828
PstI	411
SacI	834
XhoI	420

sequence landmarks (bp):

ampicillin-esistence:	2876–2016
rep (pMB1):	1856–1242
lac-promotor:	936– 897
evoglow-Bs2-gene:	830– 417

```

      EcoRI   SacI   NdeI
      |       |       |
5'  >> tac gaa ttc gag ctc cat |atg...   evoglow-Bs2   ...gag gct gca ggc atg caa gct tgg...>> 3'
      >> atg ctt aag ctc gag gta |tac...   ...ctc   ...ctc   gga cgt ccg tac gtt cga acc...>>
    
```

In pGLOW-Bs2-stop the evoglow-Bs2-gene includes a termination codon indicated in red

```

      EcoRI   SacI   NdeI
      |       |       |
5'  >> tac gaa ttc gag ctc cat |atg...   evoglow-Bs2   ...gag tga ctg cag gca tgc aag ctt ggc...>> 3'
      >> atg ctt aag ctc gag gta |tac...   ...ctc act   gac gtc cgt acg ttc gaa ccg...>>
    
```


Limitation of use:

The plasmids distributed by company evocatal GmbH under the trade mark "evoglow®" are subject to intellectual property right applications. Company evocatal GmbH commits the plasmids and genes to the customer solely for the purpose of own scientific, non-commercial investigation or own, non-commercial research and development purposes. Any commercial use requires the conclusion of a separate license contract. For commercial customers (i.e. companies or any for-profit-organisations) the permission of use for the above said research purposes of the evoglow genes and plasmids is limited to six (6) months, starting from the date of purchase. For detailed license information please contact evocatal GmbH, Germany. Furthermore, the modification of the evoglow-gene-sequences and of the plasmids is assertive prohibited. By the purchase of the plasmids without notice of defect within a term of 14 days beginning with the day of delivery, the customer agrees to the limitation of use above.

Exclusion of liability:

Company evocatal GmbH commits the plasmids distributed under the trade mark "evoglow®" to the customer solely for use in microorganism. The plasmids are not intended for the use in human medical or veterinarian medical diagnostic purposes, or in the animal or human body; any liability of the company evocatal GmbH is excluded. The committed plasmids are classified to security level S1. The national safety regulations for handling genetic engineered microorganism have to be considered and fulfilled when handling the plasmids. By the purchase of the plasmids without notice of defect within a term of 14 days beginning with the day of delivery, the customer agrees to the above exclusion of liability and indemnifies company evocatal GmbH from any liability for the use of the plasmids beyond the intended application.

Fluorescence Spectra:



Technical Data:

excitation max.:	449 nm
emission max.:	495 nm (λ_{ex} : 450 nm)
molar extinction coefficient (FMN):	$\epsilon = 12,500 \pm 500 \text{ M}^{-1} \text{ cm}^{-1}$ (at 450 nm)
quantum yield:	$Q_F = 0,39$

Sequence Information:

evoglow-Bs2 nucleotide sequence (411 bp)

```

1  ATGGCGTCGT  TCCAGTCGTT  CGGCATCCCG
31  GGCCAGCTGG  AAGTCATCAA  GAAGGCGCTG
61  GATCACGTGC  GCGTCGGCGT  GGTCATCACC
91  GATCCGCGCG  TGGAAGATAA  CCCGATCGTC
121 TACGTGAACC  AGGGCTTCGT  GCAGATGACC
151 GGCTACGAGA  CCGAGGAAAT  CCTGGGCAAG
181 AACGCGCGCT  TCCTCCAGGG  GAAGCACACC
211 GATCCGGCGG  AAGTGGACAA  CATCCGCACC
241 GCGCTGCAAA  ATAAAGAACC  GGTCACCGTG
271 CAGATCCAGA  ACTACAAGAA  GGACGGCACG
301 ATGTTCTGGA  ACGAACTGAA  CATCGATCCG
331 ATGGAAATCG  AGGATAAGAC  GTATTTTCGTC
361 GGCATCCAGA  ACGACATCAC  CAAGCAGAAG
391 GAATATGAAA  AGCTGCTCGA  G TGA

```

termination codon in pGLOW-Bs2-stop only

evoglow-Bs2 amino acid sequence (137 aa)

```

1  MASFQSFQIP  GQLEVIKKAL  DHVRVGVVIT
31  DPALEDNPIV  YVNQGFVQMT  GYETEEILGK
61  NARFLQ GKHT  DPAEVDNIRT  ALQNKEPVTV
91  QIQNYKKDGT  MFWNELNIDP  MEIEDKTYFV
121 GIQNDITKQK  EYEKLLLE

```


Note:

The nucleotide sequence of **evoglow-Bs2-stop** is identical to that of **evoglow-Bs2** but includes a termination-codon at the 3'-end as indicated in red in the adjacent sequence. The plasmid base pair numbering is hence shifted 3 positions in the 'stop'-constructs.

All plasmids in this kit are derived from pUC18 and replicate only in *E. coli* and closely related microorganisms.

Plasmids are not suitable for expression.

pGLOW-Pp1 and pGLOW-Pp1-stop
Plasmid Information

Plasmid Map


```

      EcoRI   SacI   NdeI
      |-----|-----|
5' >> tac gaa ttc gag ctc cat atg... evoglow-Pp1 ...cac gct gca ggc atg caa gct tgg c...>> 3'
   >> atg ctt aag ctc gag gta tac... ..gtg cga cgt ccg tac gtt cga acc g...>>
    
```

In pGLOW-Bs2-stop the evoglow-Bs2-gene includes a termination codon indicated in red

```

      EcoRI   SacI   NdeI
      |-----|-----|
5' >> tac gaa ttc gag ctc cat atg... evoglow-Pp1 ...cac tga ctg cag gca tgc aag ctt ggc...>> 3'
   >> atg ctt aag ctc gag gta tac... ..gtg act gac gtc cgt acg ttc gaa ccg...>>
    
```


Limitation of use:

The plasmids distributed by company evocatal GmbH under the trade mark "evoglow®" are subject to intellectual property right applications. Company evocatal GmbH commits the plasmids and genes to the customer solely for the purpose of own scientific, non-commercial investigation or own, non-commercial research and development purposes. Any commercial use requires the conclusion of a separate license contract. For commercial customers (i.e. companies or any for-profit-organisations) the permission of use for the above said research purposes of the evoglow genes and plasmids is limited to six (6) months, starting from the date of purchase. For detailed license information please contact evocatal GmbH, Germany. Furthermore, the modification of the evoglow-gene-sequences and of the plasmids is assertive prohibited. By the purchase of the plasmids without notice of defect within a term of 14 days beginning with the day of delivery, the customer agrees to the limitation of use above.

Exclusion of liability:

Company evocatal GmbH commits the plasmids distributed under the trade mark "evoglow®" to the customer solely for use in microorganism. The plasmids are not intended for the use in human medical or veterinarian medical diagnostic purposes, or in the animal or human body; any liability of the company evocatal GmbH is excluded. The committed plasmids are classified to security level S1. The national safety regulations for handling genetic engineered microorganism have to be considered and fulfilled when handling the plasmids. By the purchase of the plasmids without notice of defect within a term of 14 days beginning with the day of delivery, the customer agrees to the above exclusion of liability and indemnities company evocatal GmbH from any liability for the use of the plasmids beyond the intended application.

Fluorescence Spectra:



Technical Data:

excitation max.:	449 nm
emission max.:	495 nm (λ_{ex} : 450 nm)
molar extinction coefficient (FMN):	$\epsilon = 12,500 \pm 500 \text{ M}^{-1} \text{ cm}^{-1}$ (at 450 nm)
quantum yield:	$Q_f = 0,17$

Sequence Information:

evoglow-Pp1 nucleotide sequence (444 bp)

```

1  ATGATCAACG  CAAAACCTCCT  GCAACTGATG
31  GTCGAACATT  CCAACGATGG  CATCGTTGTC
61  GCCGAGCAGG  AAGGCAATGA  GAGCATCCTT
91  ATCTACGTCA  ACCCGGCCTT  CGAGCGCCTG
121 ACCGGCTACT  GCGCCGACGA  TATTCTCTAT
151 CAGGACGCC  GTTTTCTTCA  GGGCGAGGAT
181 CACGACCAGC  CGGGCATCGC  AATTATCCGC
211 GAGGCGATCC  GCGAAGGCCG  CCCCTGCTGC
241 CAGGTGCTGC  GCAACTACCG  CAAAGACGGC
271 AGCCTGTTCT  GGAACGAGTT  GTCCATCACA
301 CCGGTGCACA  ACGAGGCGGA  CCAGCTGACC
331 TACTACATCG  GCATCCAGCG  CGATGTCACA
361 GCGCAAGTAT  TCGCCGAGGA  AAGGGTTCGC
391 GAGCTGGAGG  CTGAAGTGGC  GGAAGTGCGC
421 CGGCAGCAGG  GCCAGGCCAA  GCAC  TGA
  
```

termination codon in pGLOW-Pp1-stop only

evoglow-Pp1 amino acid sequence (148 aa)

```

1  MINAKLLQLM  VEHSNDGIVV  AEQEGNESIL
31  IYVNPAFERL  TGYCADDILY  QDARFLQGED
61  HDQPGIAIIR  EAIREGRPCC  QVLRNYRKDG
91  SLFWNELSIT  PVHNEADQLT  YYIGIQRDVT
121 AQVFAEERVR  ELEAEVAELR  RQQGQAKH
  
```

Note:

The nucleotide sequence of **evoglow-Pp1-stop** is identical to that of **evoglow-Pp1** but includes a termination-codon at the 3'-end as indicated in red in the adjacent sequence. The plasmid base pair numbering is hence shifted 3 positions in the 'stop'-constructs.

All plasmids in this kit are derived from pUC18 and replicate only in *E. coli* and closely related microorganisms.

Plasmids are not suitable for expression.

Contact Information:

For any further information, questions or remarks please contact us:



**evocatal GmbH
Merowingerplatz 1a
40225 Düsseldorf
www.evocatal.com
info@evocatal.com**